



Comprehensive genomic profiling

Maximize the chances of identifying
molecularly matched therapies.

One biopsy, one test, one report can
lead to improved patient outcomes.



Enable precision medicine with comprehensive genomic profiling

Driven by targeted molecular therapies and immunotherapies, precision medicine offers an individualized approach that battles cancer at its core—the genome. As new oncogenic drivers are uncovered at an unprecedented rate, a testing method that can keep pace is needed. One method meeting this challenge is comprehensive genomic profiling, or CGP.



One test can make the difference

- Analyze hundreds of clinically relevant biomarkers simultaneously
- Assess DNA and RNA alterations, including SNVs, CNVs, indels, fusions, and splice variants*
- Measure TMB, MSI, and HRD^{1-5*}



Maximize data from one biopsy

- Replace multiple single-gene tests or small hotspot panels with one comprehensive test
- Decrease the need to rebiopsy to obtain more data⁶⁻⁸
- Reanalyze data as new biomarkers are discovered



Get results faster

- Achieve a faster turnaround time with CGP than sequential iterative testing⁹



Receive actionable results

- Unlock potential opportunities for molecularly matched therapy regimens
- Identify potential eligibility for matched clinical trials



Using CGP to match patients with targeted or immunotherapies has been linked to improved clinical outcomes¹⁰⁻¹⁵

- Increased ORR*
- Increased OS*
- Increased PFS*

“Comprehensive genomic profiling not only maximizes tissue utilization from often scant tumor biopsies, but it can produce biomarker information across all the biomarker classes and mutations that are required to get a comprehensive picture of the tumor...having a detailed picture of the tumor is certainly what makes decision treatment options easier for clinicians.”

Jeff Conroy
Chief Scientific Officer,
OmniSeq, and Director
of Genomics Consortium
Technologies, Roswell Park
Comprehensive Cancer Center

*SNV, single nucleotide variation; CNV, copy number variation; TMB, tumor mutational burden; MSI, microsatellite instability; HRD, homologous recombination deficiency; ORR, objective response rate; OS, overall survival; PFS, progression-free survival.

A single biomarker can transform lives

Every year, promising cancer therapies are developed and approved, offering the potential for improved clinical outcomes for patients with qualifying cancer types. Many of these therapies are molecularly targeted, zeroing in on oncogenic drivers to fight the cancer at its source. As the number of biomarkers increases, it is critical to find ways to maximize the ability to match patients with appropriate molecular treatment regimens.


Many biomarker tests are available to aid with therapy selection. However, ordering individual, sequential biomarker tests requires a significant amount of biopsy sample that is not always accessible.⁶⁻⁸ These single-gene tests screen for limited content and may miss the opportunity to identify a positive biomarker.

Identifying biomarkers is crucial

Over recent decades the underlying genomic drivers of many cancers have been revealed, creating a more complex picture of the disease than originally imagined. We now know that genetic variants can vary widely between histologically distinct tumors, and that heterogeneity can exist within a tumor itself. Moreover, tumors with similar driver mutations can have differing responses to therapies, adding additional complexity. This diversity makes it essential to understand the genetic makeup of each tumor.

COMPREHENSIVE GENOMIC PROFILING (CGP): a single assay that uses next-generation sequencing (NGS) to assess relevant cancer biomarkers, as established in guidelines and clinical trials, and inform therapy guidance. CGP detects biomarkers at nucleotide level-resolution and may comprise all major genomic variant classes (SNVs, indels, CNVs, fusions, splice variants), and large genomic signatures (TMB, MSI, HRD), maximizing the ability to find clinically actionable alterations.

Biomarkers in US guidelines and drug labels for highly prevalent tumors^{16,17}



Based on evidence in scientific literature, presence in clinical trials, or linked to labels in other histologies.
*In ovarian cancer, HRD is noted as a molecular signature based on a measurement of genomic instability.

Biomarker-driven therapy selection

Currently, 650+ drugs are listed on the National Cancer Institute cancer treatment website.¹⁸ Clinicaltrials.gov lists over 4200 biomarker-linked trials in progress, globally, for all cancer types.* With 64% of available biomarker-driven therapies for non-small cell lung cancer (NSCLC) gaining US FDA approval in just the past five years,¹⁶ it is clear that the pace of discovery is growing exponentially. The ability to detect these new biomarkers can match patients with targeted therapies, potentially leading to improved outcomes.¹⁰⁻¹⁵

Approved biomarker-driven therapies available for NSCLC treatment¹⁶

	Targeted therapies								Immunotherapy + chemotherapy
	KRAS G12c	NTRK	BRAF	ROS-1	ALK	EGFR	MET	RET	
2003						Gefitinib 3L			
2004						Erlotinib 2L			
2005									
//									
2010					Crizotinib				
2011									
2012						Erlotinib 1L			
2013					Crizotinib, 1L	Afatinib 1L			
2014						Gefitinib, 1L			Pembrolizumab, 2L (PD-L1)
2015						Necitumumab			Pembrolizumab, 1L (PD-L1)
2016				Crizotinib		Osimertinib, 2L			
2017			Dabrafenib and Trametinib		Brigatinib				
2018		Larotrectinib			Ceritinib, 1L	Osimertinib, 1L			Atezolizumab, 1L Without EGFR/ALK
2019					Alectinib, 1L	Dacomitinib, 1L			Pembrolizumab, 1L (PD-L1)
2020						Ramucirumab/Erlotinib			
2020					Brigatinib	Osimertinib, 1L	Capmatinib	Nivolumab/Ipilimumab (PD-L1), 1L	Atezolizumab, 1L (PD-L1)
2021	Sotorasib				Lorlatinib	Amivantamab-vmjw	Selpercatinib		Atezolizumab (adjuvant)
2021						Mobocertini	Tepotinib	Pralsetinib	
2022		Entrectinib		Entrectinib					

US FDA-approved indications of NSCLC treatments since 2003. Abbreviations: 1L, first-line; 2L, second-line; ALK, anaplastic lymphoma kinase; BRAF, murine sarcoma viral oncogene homolog B; EGFR, epidermal growth factor receptor; FDA, Food and Drug Administration; NSCLC, non-small cell lung cancer; NTRK, neurotropic tropomyosin receptor kinase; PD-L1, programmed-death ligand 1; ROS-1, c-ros1 oncogene.


* Based on a search at clinicaltrials.gov using a combination of the terms "genetic," "genomic," "DNA," and "RNA." Accessed May 18, 2023.

TMB and MSI predict response to immunotherapy

Clinical trials and regulatory approvals have established several cancer immunotherapy treatments for multiple tumor types.¹⁹ The ability to identify genomic signatures that help predict response to these treatments is important for better predicting patients who will benefit from these treatments. Tumor mutational burden (TMB) and microsatellite instability (MSI) status are two biomarkers that may predict patient response to immunotherapy and recommended for testing by select guidelines.



TMB measures the number of mutations in the coding sequence of the genome. In fact, 13-26% of advanced cancer patients exhibit TMB-High results across tumor types.²⁰⁻²⁴ Accurately measuring TMB requires a large panel²⁵ and cannot be performed using small targeted panels that lack sufficient gene coverage.



Not all NGS panels are the same:

Large assays with ~1.1 Mb of coding genome are needed to assess TMB accurately^{1,25}



MSI measures changes in DNA base repeats that may occur within the tumor. MSI status is FDA-approved for selection of solid tumors for treatment with checkpoint inhibitors.^{26,27} Using CGP allows for analysis of a greater number of microsatellite loci than other methods, presenting opportunities to identify new MSI profiles in previously uncharacterized cancer types.²⁸

"tbNGS [tumor-based next-generation sequencing] can provide important clinical information and has the potential to improve patient outcomes when results are effectively integrated into treatment planning."

Foster KI, et al.²⁹

About HRD

As researchers learn more about the underlying genomics of cancer, they are uncovering broader genomic signatures that may occur across cancer types. Homologous recombination deficiency (HRD) is a genomic signature of increasing importance in tumor biology for ovarian, breast, pancreatic, and prostate cancers. Patients with a tumor exhibiting HRD (HRD+) may be eligible for targeted PARPi therapies.^{30,31}

Poly (ADP-ribose) polymerase inhibitors (PARPi) are a class of targeted therapies that may be prescribed to ovarian cancer patients with tumors that exhibit HRD.³²⁻³⁵

Genomic scarring and GIS

HRD results from a cell's inability to repair double-stranded DNA breaks using the homologous recombination repair (HRR) pathway. Double-stranded breaks are either not repaired or repaired using the error-prone nonhomologous end joining (NHEJ) pathway. These alternatives can lead to genomic instability, in the form of genomic scars, resulting in tumorigenesis.³⁶

Genomic scars are aberrations that result in structural changes to the chromosomes, including loss of heterozygosity (LOH),³⁷ telomeric-allelic imbalance (TAI),³⁸ and large-scale state transitions (LST).³⁹ Each genomic scar can be measured alone, or summed together to produce a genomic instability score (GIS).

Determining HRD status

HRD status can be determined by evaluating the presence of causal genes (*BRCA* and other HRR genes) and/or the effect of genomic scarring. There is increasing evidence that assessing all three genomic scars (LOH, TAI, LST), along with causal genes, may maximize identification of HRD+ samples.⁴⁰⁻⁴²

HRD + CGP may provide added context to a disease

HRD status may reflect tumors positive for *BRCA* status and/or with a high GIS but does not address genomic variants beyond these factors. CGP can potentially identify additional biomarkers that may match patients with approved therapies or open clinical trials. By combining HRD assessment and CGP, clinicians can understand the full biology within a tumor, possibly leading to improved patient outcomes.

 To learn more, read *Powerful insights from combining HRD and CGP*⁴³

Loss of heterozygosity (LOH)



One of the two alleles for a gene is lost, creating a homozygous cell. This may result in malignant cell growth if the remaining allele does not function properly.

Telomeric-allelic imbalance (TAI)



The allele ratio at the end of the chromosomes (telomeres) in a pair do not match. That is, one chromosome has a greater number of alleles than the other.

Large-scale state transitions (LST)



Breakpoints between regions of the chromosome resulting in discrepancies within the chromosome pair.

Current testing methods may leave patients behind

Single-gene and small hotspot panel testing methods are limited in their ability to detect known and emerging biomarkers and molecular signatures, potentially missing important actionable variants.^{6,44-46} CGP provides broad coverage of the genome, capturing a comprehensive set of clinically relevant genes in one test. Based on NGS, CGP detects DNA and RNA variants, including key genomic signatures, maximizing the ability to detect actionable variants compared to conventional methods.^{6,10,25,45-49} The evaluation of RNA fusions and splice variants through RNA-based NGS assays is increasingly being recommended by clinical guidelines and associations.^{50,51}

“Single analyte tests or hotspot panels that are limited to a single gene or a few hot spot regions have really no potential to expand and keep up with the emerging markers, nor do they have the ability with their small footprint to do the larger signatures, such as MSI, TMB, and HRD.”

Jeff Conroy
Chief Scientific Officer, OmniSeq and Director of Genomics
Consortium Technologies, Roswell Park Comprehensive
Cancer Center

Small hotspot
panels can miss

81%

of actionable
biomarkers in patients
with refractory cancers,
based on study with
10,000 patients¹⁰

CGP tests can increase the number of relevant biomarkers identified compared to conventional testing approaches, such as single-gene tests and hotspot NGS panels^{6,10,45,47,48}

Single gene		EGFR	KRAS	ALK	MET	ROS	ERBB2	NTRK	RET	Other Approved	Other Emerging	TMB	MSI	HRD
	Small variants													
	CNVs													
	Fusions													
	Splice variant													

Hotspot panel		EGFR	KRAS	ALK	MET	ROS	ERBB2	NTRK	RET	Other Approved	Other Emerging	TMB	MSI	HRD
	Small variants													
	CNVs													
	Fusions													
	Splice variant													

CGP		EGFR	KRAS	ALK	MET	ROS	ERBB2	NTRK	RET	Other Approved	Other Emerging	TMB	MSI	HRD
	Small variants													
	CNVs													
	Fusions													
	Splice variant													

Comprehensive biomarker coverage to maximize ability to detect actionable variants



A missed biomarker is a potential missed opportunity

CGP provides a single test that uses minimal biopsy samples for deep analysis of biomarkers and molecular signatures linked to therapies, guidelines, and clinical trials. Data from CGP tests can be reanalyzed as new discoveries emerge, without the need to rebiopsy or to rerun the test.

With CGP, every discovery is a potential opportunity.

Multiple studies have demonstrated the ability of CGP to identify clinically relevant genomic alterations, across different tumor types

Author	Study details	Patient cohort	Percent of patient samples with actionable variants identified ^a
Wheler JJ, et al 2016 ⁵²	Prospective, single-center	339 patients profiled; refractory cancers, multiple types: ovarian (18%), breast (16%), sarcoma (13%), renal (7%), and others	93.5%
Hirschfield KM, et al 2016 ⁵³	Prospective	100 patients; diverse histology rare, or poor-prognosis cancers	94.5%
Zehir A, et al 2017 ¹⁰	Prospective	10,000 patients; advanced cancer across multiple solid tumor types	36.7%
Reitsma M, et al 2019 ⁴⁷	Retrospective	96 patients; multiple tumor types	90%
Foster K, et al 2022 ²⁹	Retrospective	409 patients; high-grade epithelial ovarian carcinoma	74.6%

a. The percent of actionable alterations identified in each study varies according to patient cohort, study type, CGP panel used, and criteria for categorizing a genomic alteration as actionable.

Consolidate testing for more insights


CGP replaces the need to run multiple independent tests with one consolidated test. This optimizes biopsy usage compared to iterative single-gene testing approaches that lead to tissue depletion and repeat biopsies.⁶⁻⁸

Identify more genomic variants from less sample

“The change from single-gene or single-biomarker testing to a comprehensive panel-based approach has been driven by a combination of factors that include inherent efficiency of a single comprehensive panel, which is key among cancers and other samples that have limited tissue.”

Dr. Jeremy Wallentine
Staff Hematopathologist and Medical Director of Molecular Testing, Intermountain Healthcare

CGP provides actionable information for therapy selection from one test, one workflow, and one report



A comparison of potential workflows for patient sample analysis without (top) or with (bottom) CGP. Without CGP testing (top), the sample is spread across multiple tests, each one yielding a separate report, or none at all for tests that may require additional sample (QNS, quantity not sufficient). The CGP workflow (bottom) requires one test with as few as five slides and generates a single report that provides information on hundreds of biomarkers and includes guidance for possible therapies and clinical trials.



Stay current with guidelines

Cancer diagnosis and treatment guidelines are updated as new biomarkers are discovered, new therapies approved, and new clinical studies released. CGP can help with the challenge of keeping abreast of these developments. Large CGP panels provide significant coverage of known and emerging biomarkers in key guidelines for multiple tumor types and genes involved in clinical trials, providing access to comprehensive data in a single analysis. In fact, large panel NGS tests are increasingly being included in clinical guidelines and recommendations for biomarker profiling in multiple solid tumor types.⁵⁴ NGS data are stored in a digital format that can be readily reanalyzed as new discoveries are made and approved, without the need for obtaining additional sample or running additional tests.


“It is highly recommended that clinical research centres perform multigene sequencing in the context of molecular screening programmes in order to increase access to innovative drugs and to speed up clinical research. This is particularly relevant in **breast, pancreatic, and hepatocellular cancers** where level II–IV alterations are numerous.”

– Mosele F, et al⁵⁴

NGS tests are increasingly recommended by clinical guidelines for biomarker profiling in multiple solid tumor types.⁵⁵ Tumor types in bold indicate TMB testing is recommended by clinical guidelines

Cancer types with NGS test recommendations		
Bladder	Gastric	Prostate
Bone	GIST	Rectal
Breast	Head and neck	Salivary
Cervical	Hepatobiliary	Small bowel adenocarcinoma
CNS	Lung adenocarcinoma	Soft tissue sarcomas
Cholangiocarcinoma	Neuroendocrine and adrenal tumors	Testicular
Colon	NSCLC	Thyroid
Colorectal	Occult CUP	Uterine
Cutaneous melanoma	Ovarian	Vulvar
Esophageal/esophagogastric junction	Pancreatic	

CNS, central nervous system; GIST, gastrointestinal stromal tumors; NSCLC, non-small cell lung cancer; CUP, cancer of unknown primary.



Not all NGS panels are the same:
Large assays with ~1.1 Mb of coding genome are needed to assess TMB accurately^{1,25}



CGP analyzes hundreds of biomarkers simultaneously, including MSI, TMB, and HRD

Guidelines recommend testing for up to

12

genomic biomarkers for NSCLC

Comprehensive genomic insights from one consolidated report


A CGP report may include a list of biomarkers identified by the tumor profiling assay and relevant information consolidated in one document. The user-curated report may include information on therapies associated with genetic variants identified through comprehensive tumor profiling. It also may indicate clinical trials linked to the detected biomarkers.^{10,11,47,56}

A single CGP report streamlines the workflow by including findings on genomic alterations specific to relevant biomarkers, targeted and immunotherapies specific to a tumor's genotype, and open clinical trials.








Therapy selection, sooner

By transforming the test-results-repeat paradigm experienced with sequential testing into a single consolidated test, CGP delivers results faster. One study describes a patient journey in which the patient went from CGP testing to therapy selection in just 8 days.⁹



The right therapy for the right patient

Molecular matching of a patient’s cancer with an approved or investigational therapy is essential for increasing the chances of achieving better clinical outcomes. Multiple studies have shown that genomically matching patients to targeted therapies or immunotherapies results in improved clinical outcomes.¹⁰⁻¹⁵


-  Increased overall response
-  Increased disease control
-  Extended progression-free survival (PFS)
-  Longer overall survival (OS)
-  Reduced treatment-related adverse effects (AE)

“Use of NGS panels in clinical practice may help to choose the best therapeutic options for the patients with actionable alterations.”

Dr. Bernard Doger de Spéville
Medical Oncologist,
Fundación Jiménez Díaz

CGP enables more patients to be matched to approved or investigational therapies.^{6,10,45,47,48} Additionally, CGP can improve patient enrollment in genomically matched clinical trials.^{10,47}

Patients are assigned a molecular therapy or eligibility for a clinical trial based on the tumor molecular profile






Matched therapy can lead to better patient outcomes

Study and description	Outcome	Reference
A retrospective review of 814 patients with NSCLC	OS for patients receiving molecularly targeted therapy: 31.8 months; 95% CI OS for patients receiving chemotherapy: 12.7 months; 95% CI	Gutierrez ME, et al ¹²
A study of ~1700 patients with advanced NSCLC	OS for patients receiving molecularly matched therapy: 18.6 months; 95% CI OS for patients receiving nontargeted therapy: 11.4 months; 95% CI	Singal G, et al ¹³
A study of 101 patients with lung adenocarcinoma who had CGP performed	Influenced treatment decision in ~50% of cases ORR: 65%	Rozenblum AB, et al ¹⁵
A study of 429 cancer patients; 62% received matched therapies	Longer PFS: Hazard ratio (HR) = 0.63; 95% CI, 0.50–0.80; P < 0.001 Longer OS: HR = 0.67; 95% CI, 0.50–0.90; P = 0.007 Higher stable disease: ≥ 6 months/partial/complete remission rate (52.1% vs 30.4% P < 0.001 nonmatched therapy)	Kato S, et al ¹⁴
A study of 149 patients with stage IV metastatic cancers; 49% received matched therapies	Longer PFS: median 6.5 vs. 3.1 months; P = 0.001; HR = 0.34; 95% CI, 0.19–0.62 Longer OS: HR = 0.42; 95% CI, 0.18–0.95	Sicklick JK, et al ⁶³

CGP can have a positive impact on patient enrollment in genomically matched clinical trials^{10,47}




11%

Rate of patients eligible to enroll in a clinical trial based on CGP test results, according to prospective study with 10,000 patients.¹⁰

Adding CGP from liquid biopsy further enables precision medicine

Circulating tumor DNA (ctDNA) in the blood can act as source material for CGP tests⁶⁴ and provide results that inform subsequent therapy selection^{65,66}, including immunotherapy.^{67,68}

ctDNA is obtained using a simple blood draw



Liquid biopsy is recommended today by guidelines

Different tumor types may have specific use cases for liquid biopsy, such as when:

- Patient is medically unfit for a tissue biopsy
- Insufficient material is available (QNS)
- Tissue biopsy is unavailable

Immuno-oncology biomarkers

CGP tests using liquid biopsy can detect genetic signatures, such as MSI and TMB:

MSI measured by CGP from blood (ctDNA) has been observed to be highly concordant with tissue-based standard-of-care testing with similar clinical utility for advanced gastric cancer patients on immunotherapy treatment.⁶⁷



TMB measured from blood (bTMB) using CGP has been associated with improved clinical outcomes when ≥ 20 mutations per megabase were detected in mNSCLC.⁶⁸⁻⁷²



Liquid biopsy is complementary to tissue

+15%

more clinically relevant mutations identified in mNSCLC when analysis from liquid biopsy is added to tissue⁷³⁻⁷⁵

Liquid biopsy for insufficient tissue quality or quantity

~20%

of the time clinical NGS is prevented from FFPE samples in advanced cancers^{73,76,77}

Clinical trials benefit from liquid biopsy

↓3x

decrease in screening time

↑2.3x

increase in enrollment rate in advanced gastrointestinal cancer compared to CGP from tissue only⁷⁸

Patient impact

Case study: More than typical back pain

An 11-year-old girl started experiencing severe back pain that was not relieved using over-the-counter medications. Upon further exploration, doctors discovered that she had a tumor wrapped around her aorta and spine.⁷⁹ Doctors operated and removed the tumor, but the cancer returned. Chemotherapy was unsuccessful. The tumors grew and the pain continued.⁸⁰

Looking for answers, doctors ordered a CGP test that combined DNA and RNA sequencing.^{79,81} CGP identified a novel *STRN-NTRK2* fusion.⁷⁹ *NTRK* fusions are extremely rare, occurring in < 1% of solid tumors.⁸² Knowing an underlying oncogenic driver of her cancer, doctors were able to enroll the girl in a clinical trial for larotrectinib.⁸³ With treatment, the tumors have disappeared and the girl no longer experiences back pain.⁸⁰



CGP identifies prevalent and rare biomarkers in a single test, maximizing the ability to identify one that is actionable.

Case study: A new standard of care for cervical carcinosarcoma

A 58-year-old woman who had previously been treated for metastatic cervical carcinosarcoma presented at the treating clinic with pain. Cervical carcinosarcomas occur in < 1% of women with invasive cervical neoplasms,⁸⁴ and therefore have not been thoroughly characterized.⁸⁵ Her initial diagnosis was treated with a combination of radiation therapies and chemotherapies. In three years, there had been no evidence of recurrence.⁸⁵

At this new visit, cancer serum biomarkers, other than human epididymis protein 4 (HE4), appeared normal. Cryosurgery was performed, during which tumor and plasma samples were obtained. These samples were subject to TMB and MSI evaluation using NGS-based mutation profiling. The results identified > 500 somatic mutations and an extremely high TMB. Targeted therapies were not available for the identified mutations; the patient was administered cryoablation followed by pembrolizumab based on the high TMB. At the time of submission of the Zhu et al manuscript, the patient had a PFS of 11 months.⁸⁵



High TMB influenced patient treatment regimen.

Realizing the promise of precision medicine

Not all cancers harbor the same variants. To add more complexity, the same variants can drive different cancer types. Just as there is not one cause, there is not one therapy. The more that is learned about the underlying genetics of a tumor, the more it becomes clear that each patient needs an individualized approach to treatment. This is precision medicine. Providing care for each patient based on the tumor's genomic makeup.

The continuous discovery of biomarkers and rampant pace of new therapy development bring us closer to achieving this goal every day. To keep pace, appropriate tests are needed. A significant change in the way cancer is characterized is imperative. We need tests that provide sensitive, accurate results quickly using minimal sample input. Actionable results that relate to therapeutic options, including clinical trials. CGP delivers this shift. A single CGP test using minimal sample can assess multiple DNA and RNA variant types across hundreds of genes. Sophisticated analytics interpret these results and produce easy-to-interpret reports that include recommended therapies and clinical trials. Using CGP, improving patient outcomes moves from a possibility to a reality.¹⁰⁻¹⁵

CGP. One sample.
One test. One report.
More opportunities.

"In the future [we need] to move from this fragmented landscape of biomarkers to a situation where, at least in Stage 4 [cancer], we do comprehensive profiling for each and every patient."

Prof. Dr. Wilko Weichert
Pathology, Technical University
Munich

A change in paradigm



References

- Chalmers ZR, Connelly CF, Fabrizio D, et al. Analysis of 100,000 human cancer genomes reveals the landscape of tumor mutational burden. *Genome Med.* 2017;9(1):34. Published 2017 Apr 19. doi:10.1186/s13073-017-0424-2
- Pestinger V, Smith M, Sillo T, et al. Use of an Integrated Pan-Cancer Oncology Enrichment Next-Generation Sequencing Assay to Measure Tumour Mutational Burden and Detect Clinically Actionable Variants [published correction appears in *Mol Diagn Ther.* 2020 Aug;24(4):505]. *Mol Diagn Ther.* 2020;24(3):339-349. doi:10.1007/s40291-020-00462-x
- Heydt C, Rehker J, Pappesch R, et al. Analysis of tumor mutational burden: correlation of five large gene panels with whole exome sequencing. *Sci Rep.* 2020;10(1):11387. Published 2020 Jul 9. doi:10.1038/s41598-020-68394-4
- Vanderwalde A, Spetzler D, Xiao N, Gatalica Z, Marshall J. Microsatellite instability status determined by next-generation sequencing and compared with PD-L1 and tumor mutational burden in 11,348 patients [published correction appears in *Cancer Med.* 2018 Jun;7(6):2792]. *Cancer Med.* 2018;7(3):746-756. doi:10.1002/cam4.1372
- Yamamoto H, Hirasawa A. Homologous Recombination Deficiencies and Hereditary Tumors. *Int J Mol Sci.* 2021;23(1):348. Published 2021 Dec 29. doi:10.3390/ijms23010348
- Drilon A, Wang L, Arcila ME, et al. Broad, Hybrid Capture-Based Next-Generation Sequencing Identifies Actionable Genomic Alterations in Lung Adenocarcinomas Otherwise Negative for Such Alterations by Other Genomic Testing Approaches. *Clin Cancer Res.* 2015;21(16):3631-3639. doi:10.1158/1078-0432.CCR-14-2683
- Lim C, Tsao MS, Le LW, et al. Biomarker testing and time to treatment decision in patients with advanced nonsmall-cell lung cancer. *Ann Oncol.* 2015;26(7):1415-1421. doi:10.1093/annonc/mdv208
- Yu TM, Morrison C, Gold EJ, Tradonsky A, Layton AJ. Multiple Biomarker Testing Tissue Consumption and Completion Rates With Single-gene Tests and Investigational Use of OncoPrint Dx Target Test for Advanced Non-Small-cell Lung Cancer: A Single-center Analysis. *Clin Lung Cancer.* 2019;20(1):20-29.e8. doi:10.1016/j.clcc.2018.08.010
- Piening BD, Dowdell AK, Weerasinghe R, et al. Comprehensive Genomic Profiling in Patients with Advanced Cancer in a Large US Healthcare System. Poster presented at: Association for Molecular Pathology (AMP) 2020; November 16-20, 2020; virtual meeting.
- Zehir A, Benayed R, Shah RH, et al. Mutational landscape of metastatic cancer revealed from prospective clinical sequencing of 10,000 patients [published correction appears in *Nat Med.* 2017 Aug 4;23 (8):1004]. *Nat Med.* 2017;23(6):703-713. doi:10.1038/nm.4333
- Soumerai TE, Donoghue MTA, Bandlamudi C, et al. Clinical Utility of Prospective Molecular Characterization in Advanced Endometrial Cancer. *Clin Cancer Res.* 2018;24(23):5939-5947. doi:10.1158/1078-0432.CCR-18-0412
- Gutierrez ME, Choi K, Lanman RB, et al. Genomic Profiling of Advanced Non-Small Cell Lung Cancer in Community Settings: Gaps and Opportunities. *Clin Lung Cancer.* 2017;18(6):651-659. doi:10.1016/j.clcc.2017.04.004
- Singal G, Miller PG, Agarwala V, et al. Association of Patient Characteristics and Tumor Genomics With Clinical Outcomes Among Patients With Non-Small Cell Lung Cancer Using a Clinicogenomic Database [published correction appears in *JAMA.* 2020 Feb 4;323(5):480]. *JAMA.* 2019;321(14):1391-1399. doi:10.1001/jama.2019.3241
- Kato S, Kim KH, Lim HJ, et al. Real-world data from a molecular tumor board demonstrates improved outcomes with a precision N-of-One strategy. *Nat Commun.* 2020;11(1):4965. doi.org/10.1038/s41467-020-18613-3
- Rozenblum AB, Ilouze M, Dudnik E, et al. Clinical Impact of Hybrid Capture-Based Next-Generation Sequencing on Changes in Treatment Decisions in Lung Cancer. *J Thorac Oncol.* 2017;12(2):258-268. doi: 10.1016/j.jtho.2016.10.021
- US Food & Drug Administration. Hematology/Oncology (Cancer) Approval Notifications. US FDA website. [fda.gov/drugs/resources-information-approved-drugs/hematologyoncology-cancer-approvals-safety-notifications](https://www.fda.gov/drugs/resources-information-approved-drugs/hematologyoncology-cancer-approvals-safety-notifications). Updated April 20, 2023. Accessed May 18, 2023.
- Pierian. Genomic Knowledgebase. www.pieriandx.com/genomic-knowledgebase. Accessed May 18, 2023.
- National Cancer Institute. Cancer Drugs. A to Z List of Cancer Drugs. NIH-NCI website. [cancer.gov/about-cancer/treatment/drugs](https://www.cancer.gov/about-cancer/treatment/drugs). Accessed May 18, 2023.
- Emens LA, Ascierto PA, Darcy PK, et al. Cancer immunotherapy: Opportunities and challenges in the rapidly evolving clinical landscape. *Eur J Cancer.* 2017;81:116-129. doi:10.1016/j.ejca.2017.01.035
- Marabelle A, Fakih M, Lopez J, et al. Association of tumour mutational burden with outcomes in patients with advanced solid tumours treated with pembrolizumab: prospective biomarker analysis of the multicohort, open-label, phase 2 KEYNOTE-158 study. *Lancet Oncol.* 2020;21(10):1353-1365. doi.org/10.1016/S1470-2045(20)30445-9
- Cristescu R, Aurora-Garg, D, Albright, D, et al. Association Between Tumor Mutational Burden Assessed by Whole-Exome Sequencing and Outcomes of Pembrolizumab Monotherapy. Poster presented at: 2020 American Association of Cancer Research Virtual Annual Meeting II; June 22-24, 2020. LB-261.
- Hellmann MD, Ciuleanu TE, Pluzanski A, et al. Nivolumab plus Ipilimumab in Lung Cancer with a High Tumor Mutational Burden. *N Engl J Med.* 2018;378(22):2093-2104. doi:10.1056/NEJMoa1801946
- Samstein RM, Lee CH, Shoushtari AN, et al. Tumor mutational load predicts survival after immunotherapy across multiple cancer types. *Nat Genet.* 2019;51(2):202-206. doi:10.1038/s41588-018-0312-8
- Goodman AM, Kato S, Bazhenova L, et al. Tumor Mutational Burden as an Independent Predictor of Response to Immunotherapy in Diverse Cancers. *Mol Cancer Ther.* 2017;16(11):2598-2608. doi:10.1158/1535-7163.MCT-17-0386
- Buchhalter I, Rempel E, Endris V, et al. Size matters: Dissecting key parameters for panel-based tumor mutational burden analysis. *Int J Cancer.* 2019;144(4):848-858. doi:10.1002/ijc.31878
- Colle R, Cohen R, Cocherneau D, et al. Immunotherapy and patients treated for cancer with microsatellite instability. *Bull Cancer.* 2017;104(1):42-51. doi:10.1016/j.bulcan.2016.11.006

27. US Food & Drug Administration. FDA grants accelerated approval to pembrolizumab for first tissue/site agnostic indication. US FDA website. [fda.gov/drugs/informationondrugs/approveddrugs/ucm560040.htm](https://www.fda.gov/drugs/informationondrugs/approveddrugs/ucm560040.htm). Updated May 30, 2017. Accessed May 18, 2023.
28. Kautto EA, Bonneville R, Miya J, et al. Performance evaluation for rapid detection of pan-cancer microsatellite instability with MANTIS. *Oncotarget*. 2017;8(5):7452-7463. doi:10.18632/oncotarget.13918
29. Foster KI, Shaw KRM, Jin J, et al. Clinical implications of tumor-based next-generation sequencing in high-grade epithelial ovarian cancer. *Cancer*. 2023;129(11):1672-1680. doi:10.1002/cncr.34724
30. Wagener-Ryczek S, Merkelbach-Bruse S, Siemanowski J. Bio-markers for Homologous Recombination Deficiency in Cancer. *J Pers Med*. 2021;11(7):612. Published 2021 Jun 28. doi:10.3390/jpm11070612
31. Label. Highlights of prescribing information—Lynparza. https://www.accessdata.fda.gov/drugsatfda_docs/label/2020/208558s014lbl.pdf. Updated May 20, 2020. Accessed May 22, 2023.
32. Ray-Coquard I, Pautier P, Pignata S, et al. Olaparib plus Bevacizumab as First-Line Maintenance in Ovarian Cancer. *N Engl J Med*. 2019;381(25):2416-2428. doi: 10.1056/NEJMoa1911361.
33. Moore KN, Secord AA, Geller MA, et al. QUADRA: A phase 2, open-label, single-arm study to evaluate niraparib in patients (pts) with relapsed ovarian cancer (ROC) who have received ≥ 3 prior chemotherapy regimens. *J Clin Oncol*. 2018;36(15):5514. doi: 10.1200/JCO.2018.36.15_suppl.5514.
34. Swisher EM, Lin KK, Oza AM, et al. Rucaparib in relapsed, platinum-sensitive high-grade ovarian carcinoma (ARIEL2 Part 1): an international, multicentre, open-label, phase 2 trial. *Lancet Oncol*. 2017;18(1):75-87. doi: 10.1016/S1470-2045(16)30559-9.
35. Stewart MD, Vega DM, Arend RC, et al. Homologous Recombination Deficiency: Concepts, Definitions, and Assays. *Oncologist*. 2022;27(3):167-174. doi: 10.1093/oncolo/oyab053
36. O'Connor MJ. Targeting the DNA Damage Response in Cancer. *Mol Cell*. 2015;60(4):547-560. doi:10.1016/j.molcel.2015.10.040
37. Abkevich V, Timms KM, Hennessy BT, et al. Patterns of genomic loss of heterozygosity predict homologous recombination repair defects in epithelial ovarian cancer. *Br J Cancer*. 2012;107(10):1776-1782. doi:10.1038/bjc.2012.451
38. Birkbak NJ, Wang ZC, Kim JY, et al. Telomeric allelic imbalance indicates defective DNA repair and sensitivity to DNA-damaging agents [published correction appears in *Cancer Discov*. 2013 Aug;3(8):952]. *Cancer Discov*. 2012;2(4):366-375. doi:10.1158/2159-8290.CD-11-0206
39. Popova T, Manié E, Rieunier G, et al. Ploidy and large-scale genomic instability consistently identify basal-like breast carcinomas with BRCA1/2 inactivation. *Cancer Res*. 2012;72(21):5454-5462. doi:10.1158/0008-5472.CAN-12-1470
40. Timms KM, Abkevich V, Hughes E, et al. Association of BRCA1/2 defects with genomic scores predictive of DNA damage repair deficiency among breast cancer subtypes. *Breast Cancer Res*. 2014;16(6):475. Published 2014 Dec 5. doi:10.1186/s13058-014-0475-x
41. Marquard AM, Eklund AC, Joshi T, et al. Pan-cancer analysis of genomic scar signatures associated with homologous recombination deficiency suggests novel indications for existing cancer drugs. *Biomark Res*. 2015;3(9). doi:10.1186/s40364-015-0033-4
42. Timms KM, Mills GB, Perry M, Gutin A, Lanchbury J, Brown R. Comparison of genomic instability test scores used for predicting PARP activity in ovarian cancer. *J Clin Onc*. 2020;38(15):1586. Published 2020 May 25. doi: 10.1200/JCO.2020.38.15_suppl.1586
43. Illumina. Powerful insights from combining HRD and CGP ebook. <https://www.illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/hrd-cgp-guide/hrd-cgp-guide-m-gl-00914.pdf>. Published 2022. Accessed May 22, 2023.
44. Suh JH, Johnson A, Albacker L, et al. Comprehensive Genomic Profiling Facilitates Implementation of the National Comprehensive Cancer Network Guidelines for Lung Cancer Biomarker Testing and Identifies Patients Who May Benefit From Enrollment in Mechanism-Driven Clinical Trials. *Oncologist*. 2016;21(6):684-691. doi:10.1634/theoncologist.2016-0030
45. Ali SM, Hensing T, Schrock AB, et al. Comprehensive Genomic Profiling Identifies a Subset of Crizotinib-Responsive ALK-Rearranged Non-Small Cell Lung Cancer Not Detected by Fluorescence In Situ Hybridization. *Oncologist*. 2016;21(6):762-770. doi:10.1634/theoncologist.2015-0497
46. Schrock AB, Frampton GM, Herndon D, et al. Comprehensive Genomic Profiling Identifies Frequent Drug-Sensitive EGFR Exon 19 Deletions in NSCLC not Identified by Prior Molecular Testing. *Clin Cancer Res*. 2016;22(13):3281-3285. doi:10.1158/1078-0432.CCR-15-1668
47. Reitsma M, Fox J, Borre PV, et al. Effect of a Collaboration Between a Health Plan, Oncology Practice, and Comprehensive Genomic Profiling Company from the Payer Perspective. *J Manag Care Spec Pharm*. 2019;25(5):601-611. doi:10.18553/jmcp.2019.18309
48. Kopetz S, Shaw K, Lee J, et al. Use of a Targeted Exome Next-Generation Sequencing Panel Offers Therapeutic Opportunity and Clinical Benefit in a Subset of Patients With Advanced Cancers. *JCO Precis Oncol*. 2019;3:1-14. doi: 10.1200/PO.18.00213
49. Samorodnitsky E, Jewell BM, Hagopian R, et al. Evaluation of Hybridization Capture Versus Amplicon-Based Methods for Whole-Exome Sequencing. *Hum Mutat*. 2015;36(9):903-914. doi:10.1002/humu.22825
50. Lindeman NI, Cagle PT, Aisner DL, et al. Updated Molecular Testing Guideline for the Selection of Lung Cancer Patients for Treatment With Targeted Tyrosine Kinase Inhibitors: Guideline From the College of American Pathologists, the International Association for the Study of Lung Cancer, and the Association for Molecular Pathology. *Arch Pathol Lab Med*. 2018;142(3):321-346. doi:10.5858/arpa.2017-0388-CP
51. Belli C, Penault-Llorca F, Ladanyi M, et al. ESMO recommendations on the standard methods to detect RET fusions and mutations in daily practice and clinical research. *Ann Oncol*. 2021;32(3):337-350. doi:10.1016/j.annonc.2020.11.021
52. Wheler JJ, Janku F, Naing A, et al. Cancer Therapy Directed by Comprehensive Genomic Profiling: A Single Center Study. *Cancer Res*. 2016;76(13):3690-3701. doi:10.1158/0008-5472.CAN-15-3043
53. Hirshfield KM, Tolkunov D, Zhong H, et al. Clinical Actionability of Comprehensive Genomic Profiling for Management of Rare or Refractory Cancers. *Oncologist*. 2016;21(11):1315-1325. doi:10.1634/theoncologist.2016-0049

54. Mosele F, Remon J, Mateo J, et al. Recommendations for the use of next-generation sequencing (NGS) for patients with metastatic cancers: a report from the ESMO Precision Medicine Working Group. *Ann Oncol*. 2020;31(11):1491-1505. doi:10.1016/j.annonc.2020.07.014
55. Data on file. Illumina, Inc., 2023.
56. Sabatini LM, Mathews C, Ptak D, et al. Genomic Sequencing Procedure Microcosting Analysis and Health Economic Cost-Impact Analysis: A Report of the Association for Molecular Pathology. *J Mol Diagn*. 2016;18(3):319-328. doi:10.1016/j.jmoldx.2015.11.010
57. Mayo Clinic Laboratories. EGFR - Specimen: EGFR Gene, Mutation Analysis, 51 Mutation Panel, Tumor. Mayo Clinic Laboratories website. www.mayocliniclabs.com/test-catalog/overview/614665. Accessed May 18, 2023.
58. ARUP Laboratories. EGFR Mutation Detection by Pyrosequencing | ARUP Lab Test Directory. ARUP Laboratories website. aruplab.com/Tests/Pub/2002440. Accessed May 18, 2023.
59. Abbott Molecular. ALK-US-CE-Clinical-PL_R3_mw001_3060.pdf. Vysis ALK Break Apart FISH Probe Kit. Abbott Molecular website. www.molecular.abbott/us/en/products/oncology/vysis-alk-break-apart-fish-probe-kit. Accessed May 18, 2023.
60. NeoGenomics Laboratories. MET Exon 14 Deletion Analysis | NeoGenomics Laboratories. NeoGenomics Laboratories website. neogenomics.com/test-menu/met-exon-14-deletion-analysis. Accessed May 18, 2023.
61. Geisinger. Specimen collection and processing instructions for BRAF MUTATION ANALYSIS. Geisinger Medical Laboratories website. www.geisingermedicallabs.com/catalog/details.cfm?tid=1740. Accessed May 18, 2023.
62. Geisinger. Specimen collection and processing instructions for KRAS MUTATION ANALYSIS. Geisinger Medical Laboratories website. www.geisingermedicallabs.com/catalog/details.cfm?tid=1638. Accessed May 18, 2023.
63. Sicklick JK, Kato S, Okamura R, et al. Molecular profiling of cancer patients enables personalized combination therapy: the I-PREDICT study. *Nat Med*. 2019;25(5):744-750. doi:10.1038/s41591-019-0407-5
64. Bettgowda C, Sausen M, Leary RJ, et al. Detection of circulating tumor DNA in early- and late-stage human malignancies. *Sci Transl Med*. 2014;6(224):224ra24. doi:10.1126/scitranslmed.3007094
65. Wan JCM, Massie C, Garcia-Corbacho J, et al. Liquid biopsies come of age: towards implementation of circulating tumour DNA. *Nat Rev Cancer*. 2017;17(4):223-238. doi:10.1038/nrc.2017.7
66. Wongchenko M, Kim S, Saura C, et al. Circulating Tumor DNA and Biomarker Analyses From the LOTUS Randomized Trial of First-Line Ipatasertib and Paclitaxel for Metastatic Triple-Negative Breast Cancer. *JCO Precis Oncol*. 2020;4:1012-1024. doi: 10.1200/PO.19.00396
67. Willis J, Lefterova MI, Artyomenko A, et al. Validation of Microsatellite Instability Detection Using a Comprehensive Plasma-Based Genotyping Panel. *Clin Cancer Res*. 2019;25(23):7035-7045. doi:10.1158/1078-0432.CCR-19-1324
68. Socinski M, Velcheti V, Mekhail T, et al. LBA83 - Final efficacy results from B-F1RST, a prospective phase II trial evaluating blood-based tumour mutational burden (bTMB) as a predictive biomarker for atezolizumab (atezo) in 1L non-small cell lung cancer (NSCLC). *Ann Oncol*. 2019;30(5):919-920. doi: org/10.1093/annonc/mdz394.081
69. Rizvi NA, Cho BC, Reinmuth N, et al. Durvalumab With or Without Tremelimumab vs Standard Chemotherapy in First-line Treatment of Metastatic Non-Small Cell Lung Cancer: The MYSTIC Phase 3 Randomized Clinical Trial. *JAMA Oncol*. 2020;6(5):661-674. doi:10.1001/jamaoncol.2020.0237
70. Garassino MC, Gadgeel SM, Rodriguez-Abreu D, et al. Evaluation of blood TMB (bTMB) in KEYNOTE-189: Pembrolizumab (pembro) plus chemotherapy (chemo) with pemetrexed and platinum versus placebo plus chemo as first-line therapy for metastatic nonsquamous NSCLC. *J Clin Oncol*. 2020;38:15_suppl, 9521-9521. doi: 10.1158/1078-0432.CCR-20-3771
71. Gandara DR, Paul SM, Kowanzet M, et al. Blood-based tumor mutational burden as a predictor of clinical benefit in non-small-cell lung cancer patients treated with atezolizumab. *Nat Med*. 2018;24(9):1441-1448. doi:10.1038/s41591-018-0134-3
72. Si H, Kuziora M, Quinn KJ, et al. A Blood-based Assay for Assessment of Tumor Mutational Burden in First-line Metastatic NSCLC Treatment: Results from the MYSTIC Study. *Clin Cancer Res*. 2020 Dec 22. doi: 10.1158/1078-0432.CCR-20-3771. Epub ahead of print. PMID: 33355200
73. Aggarwal C, Thompson JC, Black TA, et al. Clinical Implications of Plasma-Based Genotyping With the Delivery of Personalized Therapy in Metastatic Non-Small Cell Lung Cancer. *JAMA Oncol*. 2019;5(2):173-180. doi:10.1001/jamaoncol.2018.4305
74. Leighl NB, Page RD, Raymond VM, et al. Clinical Utility of Comprehensive Cell-free DNA Analysis to Identify Genomic Biomarkers in Patients with Newly Diagnosed Metastatic Non-small Cell Lung Cancer. *Clin Cancer Res*. 2019;25(15):4691-4700. doi:10.1158/1078-0432.CCR-19-0624
75. Palmero R, Taus A, Viteri S, et al. Biomarker Discovery and Outcomes for Comprehensive Cell-Free Circulating Tumor DNA Versus Standard-of-Care Tissue Testing in Advanced Non-Small Cell Lung Cancer. *JCO Precision Oncology*. 2021;5:93-102. doi:10.1200/PO.20.00241
76. Sholl LM, Do K, Shivdasani P, et al. Institutional implementation of clinical tumor profiling on an unselected cancer population. *JCI Insight*. 2016;1(19):e87062. Published 2016 Nov 17. doi:10.1172/jci.insight.87062
77. Meric-Bernstam F, Brusco L, Shaw K, et al. Feasibility of Large-Scale Genomic Testing to Facilitate Enrollment Onto Genomically Matched Clinical Trials. *J Clin Oncol*. 2015;33(25):2753-2762. doi:10.1200/JCO.2014.60.4165
78. Nakamura Y, Taniguchi H, Ikeda M, et al. Clinical utility of circulating tumor DNA sequencing in advanced gastrointestinal cancer: SCRUM-Japan GI-SCREEN and GOZILA studies. *Nat Med*. 2020;26(12):1859-1864. doi:10.1038/s41591-020-1063-5
79. Wu LW, Pavlock T, Patterson A, et al. Durable Clinical Response to Larotrectinib in an Adolescent Patient With an Undifferentiated Sarcoma Harboring an STRN-NTRK2 Fusion. *JCO Precis Oncol*. 2018;2:1-8. doi:10.1200/PO.18.00101
80. Knight B. El Paso girl survives cancer with 'miracle' surgery at UT Southwestern Medical Center. El Paso Times. Published May 10, 2018. Updated May 14, 2018. <https://www.elpasotimes.com/story/news/2018/05/10/el-paso-girl-beats-cancer-ut-southwesternmedical-center-surgery/537977002/>. Accessed May 18, 2023.
81. Foundation Medicine. FoundationOne Heme | Foundation Medicine. Foundation Medicine website. <https://www.foundationmedicine.com/test/foundationone-heme>. Accessed May 18, 2023.

82. Cocco E, Scaltriti M, Drilon A. NTRK fusion-positive cancers and TRK inhibitor therapy. *Nat Rev Clin Oncol*. 2018;15(12):731-747. doi:10.1038/s41571-018-0113-0
83. McGinley L. FDA approves 'precision medicine' drug for different cancers with same mutation. The Washington Post. November 26, 2018. https://www.washingtonpost.com/health/2018/11/26/fda-approves-precision-medicine-drug-different-cancers-withsame-mutation/?fbclid=IwAR2iksVfMf9b-70ybc9OdWS_pa3_OG7yNgg9Mdg_ScuqQ7qJ0gM2bHzlQ_M. Accessed January 19, 2021.
84. Bansal S, Lewin SN, Burke WM, et al. Sarcoma of the cervix: natural history and outcomes. *Gynecol Oncol*. 2010;118(2):134-138. doi:10.1016/j.ygyno.2010.04.021
85. Zhu B, Liu Y, Li J, et al. Exceptional Response of Cryoablation Followed by Pembrolizumab in a Patient with Metastatic Cervical Carcinosarcoma with High Tumor Mutational Burden: A Case Report. *Oncologist*. 2020;25(1):15-18. doi:10.1634/theoncologist.2019-0739



1.800.809.4566 toll-free (US) | +1.858.202.4566 tel

techsupport@illumina.com | www.illumina.com

© 2023 Illumina, Inc. All rights reserved. All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

M-GL-00129 v4.0