

# ANALYSIS OF FATTY ACID METHYL ESTERS (FAMES) USING GAS CHROMATOGRAPHY AND THE LUMA™ MULTI-CHANNEL VACUUM ULTRAVIOLET (VUV) DETECTOR

## INTRODUCTION:

The meticulous and precise analysis of Fatty Acid Methyl Esters (FAMES) constitutes a cornerstone of quality control and product safety across a multitude of sectors, from petrochemical production to pharmaceutical manufacturing and beyond. Particularly in biodiesel production, lipidomics research, and drug synthesis, the ability to accurately discern FAMES – including co-eluting isomers – is of utmost importance.

However, conventional analytical techniques frequently fall short, struggling to separate co-eluting isomers and often delivering results shrouded in uncertainty. This lack of reliability and precision can significantly affect the quality and safety of the end products, casting doubt on the efficacy of traditional approaches and accentuating the need for a more robust solution.

Enter the LUMA detector, equipped with advanced Vacuum Ultraviolet (VUV) spectroscopy and an unparalleled 12-band selectivity. LUMA's innovative design and capabilities make it an ideal tool for tackling the complexities inherent in FAMES analysis, offering a level of selectivity that allows it to overcome the co-elution of isomers and ensure accurate, reliable results.

In this Technical Brief, we delve into the myriad applications of LUMA in the realm of FAMES analysis. From its use in enhancing biodiesel quality by accurately assessing its FAMES composition, to its potential for advancing lipidomics research by enabling the identification and quantification of lipid isomers, and its promise in pharmaceutical manufacturing for accurately determining FAMES impurities in drug substances and excipients, LUMA is poised to revolutionize the landscape of FAMES analysis.

As we explore the transformative potential of LUMA across diverse industry contexts, we will delve into the ways in which it ensures reliable, accurate results, its ability to handle co-eluting isomers, and the distinct advantages it brings to both routine and complex FAMES analyses. Join us as we discover the new industry standard for FAMES analysis: LUMA.

The FAMES analysis was carried out utilizing an Agilent 8890 Gas Chromatograph coupled to a LUMA Multi-Channel Vacuum Ultraviolet Absorbance Detector powered by OpenLab CDS. The FAMES standard mixture was obtained from Supelco (CRM47885). GC-LUMA instrument method conditions for this experiment are described in Table 1.

GC Conditions		LUMA Conditions	
Injection Volume: 1µL		Makeup Gas Pressure: 14 PSI N <sub>2</sub>	
Inlet Temperature: 250°C		System Gas Pressure: 52 psi N <sub>2</sub>	
Split Ratio: Splitless 0.75 min		Flow Cell Temperature: 290°C	
Column: FAMEWAX (30m x 0.25mm, 0.25µm)		Transfer Line Temperature: 290°C	
Carrier gas: Hydrogen @ 2.5mL/min		Acquisition Rate: 5 Hz	
Oven Program: 125°C, hold 0.75 min; 5°C/min to 240°C (0 min)			

Table 1 – Instrument conditions for analysis of FAMES.

# ANALYSIS OF FATTY ACID METHYL ESTERS (FAMES) USING GAS CHROMATOGRAPHY AND THE LUMA™ MULTI-CHANNEL VACUUM ULTRAVIOLET (VUV) DETECTOR

## RESULTS:

34 FAMES were separated in under 25 minutes, coelutions are easily handled by taking advantage of the multi-Channel selectivity that LUMA provides.




Figure 1 – FAMES Standard Chromatogram split 100:1 in Band 2 (130-143nm)




Figure 2 - Single Channel (Band2) Chromatogram of Coelutions




Figure 3 – Multi-Channel (Bands 2-6) Chromatogram of Coelutions

32	Methyl lignocerate 400 µg/mL	C24:0
33	cis-4,7,10,13,16,19-Docosahexaenoic acid methyl ester 200 µg/mL	C22:26n3
34	Methyl nervonate 200 µg/mL	C24:1n9

Figure 2 - Multi-Channel Selectivity being able to differentiate and resolve coeluting FAMES 32 and 33.

ANALYSIS OF FATTY ACID METHYL ESTERS (FAMES) USING GAS CHROMATOGRAPHY AND THE LUMA™ MULTI-CHANNEL VACUUM ULTRAVIOLET (VUV) DETECTOR




Figure 3 – 12 Point Spectra of various FAMES (Methyl Stearate, cis-9-Oleic acid methyl ester, Methyl gamma-linolenate) obtained on a 100m Rt-2560 column.




Figure 4 – 24min and 12min Chromatograms of FAMES analysis leveraging LUMA multi-channel capabilities to resolve coelutions. Oven program (top image): 125°C (0.75 min), 5°C/min to 240°C (10 min); Oven program (bottom image): 125°C (0.75 min), 15°C/min to 240°C (10 min)

# ANALYSIS OF FATTY ACID METHYL ESTERS (FAMES) USING GAS CHROMATOGRAPHY AND THE LUMA™ MULTI-CHANNEL VACUUM ULTRAVIOLET (VUV) DETECTOR




Figure 5 – FAMES analysis utilizing 100m Rt-2560 column, method translated using hydrogen carrier gas.

# ANALYSIS OF FATTY ACID METHYL ESTERS (FAMES) USING GAS CHROMATOGRAPHY AND THE LUMA™ MULTI-CHANNEL VACUUM ULTRAVIOLET (VUV) DETECTOR




Figure 6 – FAMES low concentration standard (50-150ppb concentration range.)




Figure 7 – Blank and Spiked FAMES Plasma Extract Chromatograms at low concentrations (200-400ppb)

## ANALYSIS OF FATTY ACID METHYL ESTERS (FAMES) USING GAS CHROMATOGRAPHY AND THE LUMA™ MULTI-CHANNEL VACUUM ULTRAVIOLET (VUV) DETECTOR

The results of FAMES analysis using GC-LUMA were exceptional. Across a range of sample types, from pharmaceuticals, excipients, petrochemicals, and biological samples, LUMA demonstrated high sensitivity and outstanding selectivity. Its unique ability to resolve co-eluting isomers ensured accurate identifications, while its robustness and precision resulted in consistent and reliable quantification. With LUMA, FAMES analysis is straightforward and reliable, confirming its value across diverse applications.

### **CONCLUSION:**

In conclusion, LUMA has proved to be a game-changing instrument in the analysis of Fatty Acid Methyl Esters (FAMES) across various industries such as pharmaceutical, petrochemical, and lipidomics research. Its unparalleled 12-band selectivity and ability to clearly distinguish co-eluting isomers significantly outperform traditional GC methods, offering enhanced accuracy, precision, and confidence in results.

Coupled with its superior sensitivity, LUMA makes even the most complex analyses accessible, achieving reliable identification and quantification down to trace levels. Its ease of use, including streamlined data analysis and integration with existing GC systems, further underscores LUMA's appeal, making it an invaluable tool in any laboratory's arsenal.

Whether it's assessing the quality of biodiesel by accurately determining its FAMES composition, pushing the boundaries in lipidomics research, or ensuring the highest standards in pharmaceutical manufacturing by detecting FAMES impurities in drug substances, LUMA offers tangible and transformational benefits.

Ultimately, LUMA's robust capabilities and distinct advantages underscore its potential to become the new industry standard in FAMES analysis. As we move forward, we anticipate LUMA to continually redefine the analytical landscape, delivering unmatched performance and value to laboratories worldwide.

# ANALYSIS OF FATTY ACID METHYL ESTERS (FAMES) USING GAS CHROMATOGRAPHY AND THE LUMA™ MULTI-CHANNEL VACUUM ULTRAVIOLET (VUV) DETECTOR

